

Lentivirus / Retrovirus Concentration Solution (5X)

Cat. #: **P904C** FisherSci Cat.#: **NC3242189**

Volume: 100 ml, for processing 400 ml lentivirus / retrovirus supernatant (5 ml test sample is available.)

Ship at Ambient temperature **Store at** 4°C **Shelf Life:** 12 months

Quality Control: Each lot is tested for sterility and successfully concentrating Lentriviral / Retroviral particles.

Product Description

Lentivirus Precipitation Solution 5X is a mixture of polymers optimized for the precipitation of lentiviral particles. It provides a simple, fast and highly efficient method for concentrating lentiviral particles. The protocol involves mixing your lentiviral supernatant with the Lentivirus Precipitation Solution, incubate for a short period, and spin the mixture in a standard centrifuge. It can increase the lentivirus titer by up to 100 folds as quick as in 4 hours and obtain excellent recoveries without ultracentrifugation. This solution is for research use only.

- **Efficient:** up to 100 folds titer increase
- **Rapid:** < 1-hour hands on
- **Easy:** No ultracentrifugation
- **Non-toxic:** safe for all cell lines, including ES cells

Safety Precautions: Follow the recommended NIH guidelines for all materials containing BSL-2 organisms.

Protocol

1. Transfer the media containing Lentriviral / Retroviral particles from plates to a sterile vessel and centrifuge the medium at 300 x g for 10 minutes to remove cell debris.
2. Filter the supernatant through a 0.45 µm filter.
3. Transfer filtered supernatant to a sterile vessel and add 1 volume of cold Lentivirus / Retrovirus Concentration Solution (4° C) to every 4 volumes of lentivirus / retrovirus-containing supernatant. Example: 5 ml Lentivirus / Retrovirus Concentration Solution with 20 ml viral supernatant.
4. Mix well and refrigerate for 4 hours to overnight. Lentivirus / Retrovirus-containing supernatant mixed with Lentivirus / Retrovirus Concentration Solution are stable for up to 4 days at 4°C.
5. Centrifuge mixture at 1,500 x g for 30 minutes at 4°C. After centrifugation, the Lentriviral / Retroviral particles may appear as a beige or white pellet at the bottom of the vessel.
6. Discard supernatant. Spin down residual solution by centrifugation at 1,500 x g for 5 minutes. Remove all traces of fluid by aspiration, carefully, not to disturb the precipitated Lentriviral / Retroviral particles in pellet.
7. Resuspend Lentriviral / Retroviral pellets in 1/10 to 1/100 of original volume using cold, sterile PBS or DMEM at 4°C.
8. Aliquot in cryogenic vials and store at -80°C until ready for use. Done.

Related products: (see [Virus Products](#) at 101Bio.com or search by Cat. #)

Related Products	Cat. #	Cat.# at Fisher Sci.
Lentivirus 10X Titer-Up Reagent, up to 10 folds	P906	NC1792390
Retrovirus 10X Titer-Up Reagent, up to 10 folds	P909	NC1692641
Lentivirus Packaging Kit (P902 + P904P)	P904	NC1697406
293T Transfection Reagent for Lentivirus Packaging	P902	
Lentivirus Packaging Mix : ready-to-use 3rd generation HIV-based lentiviral packaging system.	P904P	
Nano-Polymer Transfection Reagent	P901	
Water Bath / Tub Disinfection Treatment Tablet, 1 cleans for 1 week.	T20	NC3352917

Cell Immortalization Kits

SV40T Antigen Cell Immortalization Kit	SV40T
hTERT Cell Immortalization Kit	hTERT
Virus Transduction Enhancer	PV100
Virus Transduction Enhancer	PV500