

## AAV Packaging Service Requisition Form

### Customer information

Name		Institute	
E-mail		Telephone	

### Packaging service request

Cat. #	Titer	Volume	Name of the virus	Quantity Needed
S503H	> 1 x 10 <sup>11</sup> GC/mL	300 µL		
S503U	> 1 x 10 <sup>12</sup> GC/mL	300 µL		
S503C	Custom, please specify your need.			

### Vector information

	Vector name	Gene size (kb)	Volume (µL)	Concentration (µg/µL)
1				
2				
3				
4				
5				

DNA prep kit Used: \_\_\_\_\_ (Endotoxin-free required!)      Buffer: \_\_\_\_\_

We require **100 µg endotoxin-free DNA** for S503H, and **500 µg endotoxin-free DNA** for S503U. All custom virus production services should come with information on the plasmids which have to be packaged. An electric version of sequencing information (Vector NTI or genebank) is needed. Vector must be able to produce virus (no internal poly A signal, no toxic genes, no unusual 2nd structure, or 5' to 3' ITR size over 4.5kb.)

### Other requirement

Please send filled form by email ([support@101Bio.com](mailto:support@101Bio.com)) or fax (1-888-871-2390). We will response within 1 business day.